

Phase 1: Environment & Data Preparation

Before running the code, ensure your folder structure is initialized and you have "seed" data for the model to process.

1. Install Dependencies:

Open your terminal in the project root and run:

```
pip install flask torch torchvision opencv-python pillow matplotlib seaborn scikit-learn
```

2. Initialize Folders:

Ensure the following directory structure exists:

Plaintext

/static/uploads/

/static/results/

/data/labeled/

/data/unlabeled/

/models/

3. Generate Test Data:

Run your `setup_data.py` script. This creates dummy images so the training pipeline doesn't crash due to empty folders.

Terminal

```
python setup_data.py
```

Phase 2: The Training Pipeline (Algorithms 4 & 5)

This is where the **SimCLR** and **Semi-Supervised** logic happens.

1. Run the Pipeline:

Execute `python train_pipeline.py`.

2. What happens internally?

- **SimCLR Pre-training:** The model looks at images in `/data/unlabeled/` and learns to group similar retinal features.
- **Fine-tuning:** The model uses the small set in `/data/labeled/` to learn the 5 DR classes.

- **Pseudo-labeling:** The model makes high-confidence guesses on the remaining unlabeled data to improve itself.
3. **Output:** You will see a file named `dr_model_final.pth` appear in the `/models/` folder. **Do not proceed until this file exists.**

Phase 3: Analytics Generation

To populate the dashboard with the "Reports" you requested, you must generate the performance graphs.

1. **Execute Reporter:**

Terminal

```
python generate_reports.py
```

2. **Verification:** Check `static/results/`. You should now see `training_report.png`, `confusion_matrix.png`, and `tsne_clusters.png`.

Phase 4: Launching the Dashboard (Inference & XAI)

Now that the model is "smart" and the reports are ready, launch the web interface.

1. **Start Flask:**

Terminal

```
python app.py
```

2. **Access the UI:**

Open your web browser and go to: `http://127.0.0.1:5000`

Phase 5: Using the System

1. **Upload:** Select a retinal fundus image (from your data folder or a sample image).
2. **Analyze:** Click "**Start AI Analysis**".
3. **Interpret Output:**
 - **Classification:** View the DR grade (e.g., "Moderate DR").
 - **Confidence:** Check the percentage bar (e.g., "94.5%").
 - **Grad-CAM Heatmap:** Look at the colored overlay on the eye image. Red spots indicate where the model detected hemorrhages or lesions.

Summary of Script Execution Order:

Step	Script	Purpose
1	setup_data.py	Creates folders and dummy images.
2	train_pipeline.py	Executes SimCLR and Pseudo-labeling.
3	generate_reports.py	Generates the performance graphs (Analytics).
4	app.py	Starts the Explainable AI Dashboard.